



PCO-L-01 (Calza larga de Lycra)


TALLE	CIRCUNSFERENCIA DE CADERA (17)	CINTURA - TOBILLO (VI)	INGLE - TOBILLO (VII)
1	49	48	30
2	52	50	32
3	55	56	36
4	58	58	38
5	61	64	42

**Custom measurement form:
PCO ORTHOSES CARD**

PATIENT ID _____

PRODUCT CODE _____

NOTES _____




! WARNING!

- Chest and upper limb measurements should be taken with patient lying down.
- Hip and lower limb measurements are taken with patient standing.
- All measurements should be taken in centimeters!
- If the finger put between the child's body and the PCO product meets resistance, the size of the product has been selected correctly.

IMPORTANT!

Only a correctly filled in measurement card is a guarantee of proper compression of the product. We are not responsible for wrongly taken measurements!

HOW TO MEASURE?



- I** Arms circumference at the level of the head of the humerus
- II** Shoulder circumference measured through the armpit and acromion
- III** Length measured from armpit to wrist
- IV** Height measured from the armpit to the greater trochanter of the femur
- V** Measurement from the manubrium of the sternum through the perineum to the cervical vertebra C7
- VI** Height measured from waist to ankle
- VII** Height measured from parineum to ankle
- VIII** Patient is height