

CE FDA

PrimeCast™

Soft



Prime
Medical Inc.

PrimeCast SOFT es una venda plástica semi-rígida que al momento de cambiar sus propiedades inmoviliza el miembro del cuerpo en donde fue aplicado. La rigidez se basa en el números de capas utilizadas (2-3 capas dan flexibilidad, 5-6 capas rigidez). Esta venda es muy fácil de aplicar y reducen el tiempo de rehabilitación.

Orden de pedido

N°. Ref.	Descripción	Tamaño	Rollos/ Caja Interna	Caja Interna/Caja Externa
S2	PrimeCast™ Soft 2inch	5.0cm x 3.6m (2in x 4yd)	10Rolls/Box	9Boxes/CTN
S3	PrimeCast™ Soft 3inch	7.6cm x 3.6m (3in x 4yd)	10Rolls/Box	9Boxes/CTN
S4	PrimeCast™ Soft 4inch	10.2cm x 3.6m (4in x 4yd)	10Rolls/Box	9Boxes/CTN
S5	PrimeCast™ Soft 5inch	12.7cm x 3.6m (5in x 4yd)	10Rolls/Box	9Boxes/CTN

Available Colors: White, Yellow, Orange, Blue, Pink

Prime
MEDICAL INC

Prime Medical Inc. Korea
189 Gumdang-ri, Mado-myon, Hwasung-si,
Gyeonggi-do, Korea
Tel. +82-31-355-8374 / Fax. +82-31-355-8320
sales@primemedical-inc.com
www.primemedical-inc.com

Prime Medica Promotion Pvt. Ltd. India
107, Neelkanth Estate No. 1 Sativali,
Vasai (E), Dist: Thane-401 208 Maharashtra, India
Tel. +91-250-3291725 / 6526135
Fax. +91-250-2480551
india@primemedical-inc.com